



Buzzaround 4-Wheel GB152



The new Buzzaround XL+ 4-wheel scooter features compact portability with a robust 325 lb. weight capacity and comfort spring suspension for a smooth ride. It features a 3.5" ground clearance, up to 9.2-mile operating range on the standard 12 AH lithium-ion battery that weigh only 7 lbs., and a contoured stadium style mesh seat for more comfort and increased airflow. The GB152 offers a stylish design with fresh features like an EZ pivot tiller adjustment system that can be changed by hand or foot, a full-color dash with speed and battery indicators, and customizable LED side accent lighting with over 60 color configurations! The XL+ offers our new MicroTurn™ Technology with an impressive 39" turning radius to navigate even the tightest spaces with ease!



FDA Class II Medical Devices designed by Golden Technologies are intended to assist persons with chronic physical disabilities.

Standard Colors:



Velvet Red
(Installed)



Coastal Blue
(Included in Box)

Color Through Shrouds!

Customizable LED Side Accent Lighting!

Airline-Friendly Lithium-ion Battery!

MicroTurn™ Technology!

Standard Features:

- 325 lb. Weight Capacity
- Turning Radius of a 3-W with Stability of a 4-W
- 3.5" Ground Clearance
- Color LED Console with Battery Meter & Speedometer
- Front & Rear LED Lighting
- Customizable LED Side Accent Lighting with Over 60 Color Configurations
- Front & Rear Suspension
- EZ Connect Magnetic Charging Port & Cable
- Width & Angle Adjustable Flip-up Armrests
- USB A&C Ports on the Tiller
- Phone Holder
- Front Basket with Carry Handle & Drink Holder
- Quick & Easy Disassembly





LIGHTWEIGHT PORTABLE SCOOTER

Item	GB152
Description	XL ⁺
Medicare Code	K0801
Weight Capacity	325 lbs.
Drive Wheels	Rear
Maximum Speed	4.7 mph
Operating Range up to ¹	9.2/18 miles
Under Deck Clearance	3.5"
Turning Radius	39"
Type Battery	12AH/24AH Lithium-Ion
Colors	Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off-board
Controller	50 Amp
Dimensions	
Length	42"
Width	22"
Height (ground to top of back of seat)	35.5"-37.5"
Ground to Top of Seat	21.5"-23.5"
Ground to Top of Deck	6"
Deck to Top of Seat	14.5"-16.5"
Number of Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	30.6"
Weight of Unit: (assembled)²	97 lbs.
Front Body	42.5 lbs.
Rear Body	28 lbs.
Seat w/Arms	16 lbs.
Battery	7 lbs
Tires:	PU Foam
Front	8"
Rear	8"
Rear Anti-Tip Wheels	2"
Standard Seat	Stadium Mesh Back
Seat Back Height - no headrest	17"
Width x Depth	16"W x 15.5"D
Color	Black
Warranty:	
Frame	Lifetime Warranty
Drive Train	2 Year Warranty
Electronics ²	2 Year Warranty
Battery	1 Year Warranty



EZ Pivot Tiller System That Can Be Adjusted By Hand or Foot



Customizable LED Side Accent Lighting with Over 60 Color Configurations

microTURN



NEW MicroTurn Feature on GB152 Offers an Impressive 39" Turning Radius



Full-color LED Console with Battery & Speed Indicators



EZ Connect Magnetic Charging Port & Cable



Disassembles into 5 Pieces!



¹ Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain, and battery type.
² Electronics warranty excludes battery.