

Bariatric Sentra EC Heavy-Duty Wheelchair

HCPCS: K0007,K0007



- Weight Capacity: 450 lbs
- Carbon Steel frame with powder-coated silver vein, chip-proof, maintainable finish
- Double embossed vinyl upholstery is durable, lightweight, attractive and easy-to-clean
- Standard dual axle provides easy transition of seat height to hemi-level
- Polyurethane tires mounted on composite wheels provide durability and low maintenance
- 8" front casters, adjustable in three positions
- Maintenance free, sealed, precision bearings
- Padded armrests provide additional comfort
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Aluminum foot plates
- Comes with push-to-lock wheel locks

PRODUCTS	Description
STD20ECDDAHD-ELR	SentraEC,20"HD,DetDeskArm,ElvLg,1/bx
STD20ECDDAHD-SF	SentraEC,20 HD,DetDeskArm,SwgFt,1/bx
STD20ECDFAHHD-ELR	SentraEC,20in HD,DetFullArm,ElvLg,1/bx
STD20ECDFAHHD-SF	SentraEC,20"HD,DetFullArm,SwgFt,1/bx
STD22ECDDA-ELR	Sentra EC,22"Det.DeskArm,ElevLg,1/cs
STD22ECDDA-SF	Sentra EC,22in Det.DeskArm,SwgFt,1/cs
STD22ECDFAH-ELR	Sentra EC,22in Det.FullArmElevLg,1/cs
STD22ECDFAH-SF	Sentra EC,22"Det.FullArmSwgFt,1/cs
STD24ECDDA-ELR	SentraEC,24"DetDskArm,ElevLgRest,1/cs
STD24ECDDA-SF	SentraEC,24in DetDskArm,SwgFt,1/cs
STD24ECDFAH-ELR	Sentra EC,24in Det.FullArmElevLg,1/cs
STD24ECDFAH-SF	Sentra EC,24in Det.FullArmSwgFt ,1/cs

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	56 lb 57 lb 58 lb
Product Dimensions	42" (L) w/ riggings x 28" (W) x 36" (H)
Seat Dimensions	20" (W) x 18" (D) 22" (W) x 18" (D) 24" (W) x 18" (D)
Weight Capacity	450 lb
Seat To Floor Height	17.5" - 19.5"
Carton Shipping Weight	62 lb 63 lb 64 lb
Folded Dimensions	42" (L) w/ riggings x 12.5" (W) x 36" (H)

SPECIFICATIONS**Description**

Back Height	15.75"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]

ACCESSORIES**Description**

STDS818	Anti-Tippers,rear,SentraRecliner,1pr/bx, BX
STDS801	Brake Handle Ext.,8", 1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS855	Seat Belt,60" Auto-type,1/ea, EA
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
STDS834	OverheadAnti-Theft Bar,Single,1/ea, EA
STDS820	Telescoping IV Pole,W/C, 1/bx, BX
STDS803	Oxygn Tank Holdr for W/C,1/cs, CS
STDS831	Heel Loops, 1pr, BX
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
WASR	Limb Support, Right,1/cs, CS
WASL	Limb Support,Left,1/cs, CS